



*Established 1973*

Group IV Dental Lab, Inc.  
1585 SW Marlow Suite 105  
Portland Oregon 97225

503-297-2431 phone 503-297-4281 fax  
[groupivdental@gmail.com](mailto:groupivdental@gmail.com)  
[www.group4dentallab.com](http://www.group4dentallab.com)

Date \_\_\_\_\_ Deliver to Clinic Date \_\_\_\_\_

Doctor \_\_\_\_\_ License Number \_\_\_\_\_ Phone Number \_\_\_\_\_

Address \_\_\_\_\_

Patient Name \_\_\_\_\_ Gender \_\_\_\_\_ Age \_\_\_\_\_

Tooth Number \_\_\_\_\_ Shade \_\_\_\_\_ Stump Shade \_\_\_\_\_

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
32	31	30	29	28	27	26	25	24	23	22	21	20	19	18	17

### Items Included with Case

- Master Impression
- Model of Provisional
- Bite Registration
- Analog
- Other \_\_\_\_\_
- Opposing Impression or Model
- Articulator
- Photos
- Impression Coping
- Diagnostic Wax-Up
- Facebow Transfer Jig
- Photo Card
- Abutment
- Pre-Op Model
- Stick Bite
- Flash Drive
- Final Screw

### Type of Restoration Desired

- Porcelain Fused to Metal (PFM)
  - High Noble  Noble  Chrome Cobalt
- e.Max - Monolithic
- e.Max press over Zirconia
- Prettau (Monolithic Zirconia / MZ)
- Porcelain Butt Margin
- Other \_\_\_\_\_
- Full Gold Crown (FGC)
  - High Noble  Noble
- e.Max - Layered
- Porcelain Fused to Zirconia
- PMMA Temporary (polymethyl methacrylate)
- Diagnostic Wax up

### Case Notes and Goals

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

DOCTOR'S SIGNATURE \_\_\_\_\_

# Implant Restorations

## Screw Retained

- PFM  Layered Zirconia  e.Max press over Zirconia  Implant Temporary

## Cement Retained

- Custom Titanium Abutment  Custom Zirconia Abutment  
 Titanium Hue  Gold Hue  Zirconia Connection \*  Metal Connection \*\*  
 Custom Cast Gold Abutment  Other \_\_\_\_\_  
 Yellow High Noble  Yellow Noble

## Milling Center Preferred

- Group IV  Nobel Biocare  Atlantis  Straumann  3i  No Preference


## If Inadequate Clearance:

- Provide Reduction Coping  Reduce Prep without Reduction Coping  Reduce Opposing  Metal/Zirconia Island  Please Call Doctor


## Return Case for:

- Die Trim  Doctor Evaluation after Mounting  Framework Try-in  Bisque Bake

### FRAMEWORK DESIGN



### ADDITIONAL SHADE NOTES



\* Implants with zirconia connections will be manufactured by Nobel Biocare, Straumann, or Atlantis

\*\*Implants with metal connections on zirconia abutments can be milled by Group IV or Nobel Biocare.

Metal connections are cemented to the zirconia abutment in lab.

Pricing varies by manufacturer, most companies can produce platforms on their competitors, call us or consult your rep for more information. Group IV warranties zirconia custom abutments for breakage or defects in workmanship for 5 years. No warranty is provided for the implant fixture itself. Some implant companies will negate their warranty on the implant if the abutment used is manufactured by a competitor, partner, or third party. Group IV assumes no responsibility for the implant fixture. Please consult your warranty information provided by the implant manufacturer.